

Opn for Infants and Children with Hearing Loss

Dave Gordey, PhD
Director of Clinical Research &
Professional Relations, Oticon
A/S



Children live and learn in a beautifully chaotic world

Being able to make sense of sound is crucial to children's development and growth

“

*On average, children spent **80% of their total time in a mixture of speech in noise** ... and seldom were in situations classified as quiet, speech alone, or noise alone”*

Cruckley et al., 2011



Listening challenges of children with hearing loss

Learning and development



Noise impairs speech understanding and word learning

(Riley & McGregor, 2012)



Learning is negatively affected by noise and hearing loss – a **double disadvantage**

Well-fitted hearing aids have a positive influence on language development

(e.g. Pittman 2017; Tomblin et al., 2015)



In Noise, Omni Directional Technology may Overload Young Minds



Traditional Directional Technology Closes Down Sound



Oticon Opn Play and Open Sound Navigator a breakthrough in paediatric hearing care





Oticon Opn Play™

Oticon
BrainHearing™
Technology



Powered by
Velox S™



OpenSound
Navigator



Speech Guard
LX



YouMatic
LX



Speech
Rescue LX



Spatial
Sound LX



Clear
Dynamics



OpenSound
Optimizer



Wind Noise
Management







**OpenSound
Navigator**

Open Sound Navigator with 360° access to sound

- 1: Analyze
- 2: Balance
- 3: Noise removal

An aerial photograph of a park scene with various people and objects. A large, semi-transparent green circle is centered in the image, representing a 360-degree sound field. Inside this circle, there are several small, circular icons connected by lines to a central point, indicating sound sources or analysis points. These icons include a dog, a car, a person, a group of people, a stroller, and a truck. Outside the circle, there are more people and a truck. In the bottom right corner, there is a small icon of a person's head with sound waves emanating from it.

- # Open Sound Navigator with 360° access to sound
- 1: Analyze
 - 2: Balance
 - 3: Noise removal
- 
- The image shows an aerial view of a park with various people and objects. A large, semi-transparent green circle is centered in the park, representing a 360-degree sound field. Inside this circle, there are several small, circular icons connected by lines to a central point, indicating sound sources or analysis points. These icons include a dog, a car, a person, a group of people, a stroller, and a truck. Outside the circle, there are more people and objects, including a man walking a dog, a woman pushing a stroller, a group of people sitting on a picnic blanket, a man standing next to a truck, and a woman walking. The overall scene is a vibrant, sunny day in a park.



OpenSound Navigator



Analyze




Balance



Noise removal



Noise Removal





OpenSound
Optimizer

Open Sound Optimizer: Rethinking Feedback!




From FEEDBACK
MANAGEMENT TO
FEEDBACK
PREVENTION



*For prescribed fittings, according to best practice



Open Sound Optimizer

Patented
Breaker Signal



Up to 25% more speech cues to the brain

Oticon Opn Play breaks a law of physics with OpenSound Optimizer





Evidence Supporting OSN



Tween and Teen Focus Group: What do they want from their hearing devices?



OpenSound Navigator™ and Benefits for Speech Understanding in Noise for Children.

The influence of OpenSound Navigator on Listening Effort and Speech Recognition in Noise for Children.



The Usability and Listening Experiences of Tweens and Teens using the Oticon Opn and ConnectClip Microphone.

Research Studies

Building Evidence Supporting the Use of Opn for Children and Teens with Hearing Loss

- Tween and Teen Focus Group:
- What do they want from their hearing devices
- (Spangler and Gordey, 2016)

Understanding
barriers and
facilitators to
consistent
hearing device
use




If you could create the perfect hearing aid, what would it look like and what are the things it could do?

- ▶ It would have batteries that would last longer/rechargeable.
- ▶ Lightweight, slim, small.
- ▶ Eliminate all feedback.
- ▶ All controlled by my phone.




If you could create the perfect hearing aid, what would it look like and what are the things it could do?

- ▶ Better retention options for sports; sports mode for hearing aid so they can work well underneath a hockey or football helmet.
- ▶ Tracking device via an app in case they get lost.
- ▶ Hearing aid could connect to electronics, computers, internet.



**OpenSound
Navigator and
Benefits for Speech
Understanding in
Noise for Children**

**Boys Town National
Research Hospital,
Omaha NE**



Benefits of OpenSound Navigator™ (OSN) in Oticon Opn™ Hearing Aids for Children



Based on the study: "Integrated directionality and noise reduction: Effects on Children's masked thresholds"
by Jenna M. Browning, Emily Buss, Mary Flaherty, Lori J. Leibold.

Human Auditory Development Laboratory, Boys Town National Research Hospital, Omaha, NE,
and Department of Otolaryngology/Head and Neck Surgery, The University of North Carolina at Chapel Hill, Chapel Hill, NC

A realistic test environment replicating real life

Scenario 1: Facing speaker vs facing away

Purpose:

To investigate the benefit of OSN on speech understanding in noise when the children directly faced and faced away from the target speech

Participants:

- ▶ 14 participants with hearing impairment
- ▶ Mean age = 10 years
- ▶ Mild to severe sensorineural hearing loss


Stimuli in free field:

- ▶ Word recognition test
- ▶ Target speech from 0° or 60°
- ▶ Masking noise: Steady-state noise


Scenario 1:

Facing speaker vs facing away in steady-state background noise

Facing



Facing away




Results

Equal benefit regardless of facing or facing away from the talker

- ▶ OSN gives around 4 dB of benefit regardless of facing or facing away
- ▶ Children do not need to face the talker directly (unlike traditional directionality)

Up to 30% improvement in speech recognition



A realistic test environment replicating real life


Scenario 2: Multi-talker speech

Purpose:

To investigate whether OSN preserves competing speech coming from different directions behind the listeners

Stimuli in free field:

- ▶ Word recognition test
- ▶ Target speech from 0°
- ▶ Masking noise: Speech (2-talker babble)




Oticon whitepaper, 2017

Results

Speech recognition

- ▶ Speech recognition in omni-directional technology versus in OSN did not differ statistically

**OSN preserved
competing speech coming
from behind**



OpenSound Navigator offers optimal conditions to listen and learn

OpenSound Navigator offers significant benefits compared to traditional technology regarding

- Speech understanding (with or without facing the speaker directly)
- Preserving competing speech, hence allowing children to overhear

Study

1




The influence of OpenSound Navigator on Listening Effort and Speech Recognition in Noise for Children

**VU University Medical Center, Amsterdam,
The Netherlands**

Study: OpenSound Navigator – Effort

Background



Purpose:

Investigate the influence of OpenSound Navigator on listening effort and speech recognition in noise for children.

Participants:

- ▶ 10 participants with hearing impairment
- ▶ Mean age = 14.1 years old
- ▶ Mild to severe hearing loss

Oticon whitepaper, 2018 (in prep)

Study: OpenSound Navigator – Effort


Methods

Stimuli in free field:

- ▶ Sentence recognition in noise
- ▶ Listening effort
 - ▶ Self-reported (subjective measure)
 - ▶ Pupillometry (objective measure)
- ▶ Target speech from 0°
- ▶ Masking noise: Speech and steady-state noise

Sound processing:

- ▶ Prescription target:
DSL v5.0




Results

Speech recognition

- ▶ OSN improved speech recognition in noise in both listening environments by up to **5 dB SNR**

**Results are
consistent with
previous findings**



Results


Listening effort

Subjective listening effort

- ▶ Children perceived less effort while listening to speech in noise with OSN

Objective listening effort

- ▶ Pupillometry results showed a tendency that OSN slightly reduced the average pupil response in the simple listening environment, indicating less listening effort



Open sound experience with less effort

OpenSound Navigator offers improved benefits compared to traditional technology regarding

- Less effort while listening in noise
- May allocate more resources to acquire new skills and other activities

Study

2



In collaboration with
Cincinnati Children's
Medical Center



Thin Tube Hearing Aids in Tweens and Teens: User Experiences with the Oticon OPN and Connect Clip Remote Microphone

Aim(s)

- ▶ Learn more about the objective and perceived changes with OPN processing versus previous hearing processing strategies
- ▶ Comparing previous hearing device technology to receiver-in-the-ear or slim-tube technology
- ▶ Comparing the use of the Connect Clip remote microphone to traditional or Roger FM




Study Design

- Mixed Methods
- Purposeful sampling
- Survey
- Interviews




Study Design




Inclusion Criteria (n=20)

- ▶ Typical developing 12-18 year olds (high school or lower)
- ▶ Bilateral sensorineural hearing loss that fits within the OPN fitting range
- ▶ Previous use of remote microphone technology but does not have to be currently using
- ▶ No contraindication for receiver-in-the-ear technology
- ▶ Willingness to utilize Connect Clip in educational setting
- ▶ No previous use with OPN technology/processing




Preliminary Results (n=8) Survey

Q10 With my current hearing aids, my understanding of speech in noisy places is (select one):




Preliminary Results (n=8) Survey

Q11 Compared to my older hearing aids, my understanding of speech in noisy places with the new hearing aids are (select one):



Q31 Overall, how satisfied were you with the ConnectClip? How did it compare with your old FM system?




Preliminary Results (n=6)

Skype/Phone Interviews

Emerging Themes

- ▶ Inclusion
- ▶ Confidence
- ▶ Comfort
- ▶ Access
- ▶ Ease of Use



Participant Quotes

"I never knew music had so many layers! I'm hearing parts of the music I've never heard before!"

"I felt like I didn't have to really ask people around me to repeat what they said more often, which made me feel a lot better, and it just made me feel less ... I guess you could say anxious about it."

"Really, it was just being able to hear better in a noisy environment, like a crowded cafeteria or a crowded store. It was a lot easier to hear, then. Makes life easier."


"Well it makes me more comfortable with speaking to people and answering them when they ask me something. I know that I'm saying the right thing instead of just guessing because I don't really know what they said."

**12 colors to
meet the
needs of
children**



And decorative stickers to accompany them

For more personality



OTICON | **Opn Play**



Oticon Opn Play is Rechargeable



Designed to stand up to the test of childhood

Offering the best hearing care without compromising safety and wellbeing

- Robustness
- IP68 certified
 - Also with adapter drawer
- Nano-coated
- Tamper-resistant battery drawer
- LED light
- Free of allergens



A Full Range Of Easy-to-use Solutions For The Classroom



TwinLink
NFMI + 2.4 GHz

**Completely
integrated solution**




Other systems

**Universal
ear-level
solution**



**General
neckloop solution**



**Design
integrated solution**

Amigo FM

Summary

New Velox S platform is more advanced and powerful than ever


OpenSound Navigator lets children enjoy a full and balanced soundscape

OpenSound Optimizer optimizes audibility for better language development

Solid evidence package

A broad **range of connectivity** options for enhanced learning and communication





Let every child open up to the world

– even in complex learning environments



Thank you
davg@oticon.com