

Science made smarter

Callisto™

The clinic
made
portable



Redefining
portability.
- Audiometry, REM,
and Visible Speech
Mapping in one box.
Offering HIT with the
TBS10 test box.



Interacoustics

Audiometry
~~~~~

Tympanometry

ABR

OAE

Hearing Aid Fitting  
~~~~~

Balance

Set your own combination

on one user-friendly platform

Customize your Callisto™ through dedicated software modules to create your personal combination for Audiometry, REM, HIT with the TBS10 test box, and Visible Speech

“Inspired by advanced Audiometry, Real-ear Measurement & Hearing Instrument Testing systems, but wrapped in a small package”

AC440

Audiometry

REM440

Real ear measurement

VSP440

Visible speech mapping

HIT440

Hearing aid testing

Viot™

Video otoscopy

The Callisto™. Small footprint and lightweight (565g/1.25lbs)



TBS10 test box for Hearing Instrument Testing



Viot™ video otoscope offers full integration in the Callisto™ Suite

Fitting into a modern practice

The Callisto™ is a portable PC-based modular system that can include both a diagnostic audiometer, a full featured Real-Ear Measurement module (including Visible Speech Mapping) and a Hearing Instrument Testing module.

Highlights

- Intuitive - combines advanced testing and counseling
- Conforms to other Interacoustics Suite Software reducing the learning curve
- Mix and match Audiometry, HIT, REM and VSP
- Tailored customizable test routines
- Merge Video Otoscopy into the suite
- NOAH/OtoAccess™ compatible

Bring the clinic to your client

- Ideal for travelling and in-home visits
- Small footprint and lightweight (565g/1.25lbs)
- Powered via the PC USB
- Optional dedicated carrying bag on wheels
- Leave Callisto™ in the bag while testing

Audiometry


- Customizable test setups, printouts and report templates
- Visual and auditory counseling tools (Speech Banana, Hearing Loss Simulator, Master Hearing Aid)
- Digitized word lists from your PC hard-drive
- Audiogram historical comparison overlay
- High Frequency Audiometry
- Masking Help
- TEN test
- QuickSIN (optional)
- Customizable counseling overlays

Real-ear measurements


- Visible Speech Mapping with percentile analysis
- Multiple speech stimuli, in several languages, including ISTS and Live Voice
- On Top mode for hearing aid fittings
- NAL-NL2 and DSLv5
- Open fit capable
- Directionality testing
- Switch between HL and SPL views
- Coupler Base option for Coupler Based Verification
- Binaural REM
- Compatible with AutoFit, REMfit and EXPRESSfit

Hearing Instrument Testing

- Frequency Range: 100 - 10,000 Hz
- Built-in test sequences based on IEC and ANSI standards
- Test sessions can be saved and recalled later for direct comparison
- TBS10 Test Box
- Coupler based hearing aid fitting (RECD inside TBS10 Test Box)



Audiometry with counseling overlay



Real Ear Measurement (REM)



Visible Speech Mapping (optional)

Science made smarter

Interacoustics is more than state-of-the-art solutions

Our mission is clear. We want to lead the way in audiology and balance by translating complexity into clarity:

- Challenges made into clear solutions
- Knowledge made practical
- Invisible medical conditions made tangible and treatable

Our advanced technology and sophisticated solutions ease the lives of healthcare professionals.

We will continue to set the standard for an entire industry. Not for the sake of science. But for the sake of enabling professionals to provide excellent treatment for their millions of patients across the globe.

Interacoustics.com

Interacoustics A/S

Audiometer Allé 1
5500 Middelfart
Denmark

+45 6371 3555
info@interacoustics.com

interacoustics.com

Go online to
explore our
full product
range

Related products



Affinity^{2.0}
The integrated fitting
solution



TBS10
Test box



Viot™
Video otoscope

Product specifications

All technical and hardware specifications concerning all products can be downloaded from our website.



Interacoustics